

ACQREF™ INTRODUCER SHEATH



VASCULAR ACCESS


For Auxiliary EP catheter placement
(7 F)

INTEGRATED ELECTRODE

Provides a stable, indifferent
electrode that can be used for
unipolar pacing and recording

DUAL PURPOSE INTRODUCER



Eliminates the need for an additional
introducer and quadrapolar catheter



Compatible with the AcQMap® High Resolution Imaging and Mapping System

ACQREF 7F INTRODUCER FOR UNIPOLAR PACING AND RECORDING

- No inconsistent WCT
- No separate stick
- No need for separate catheter in the IVC
- Stable Electrode
- Not in contact with other catheters or the myocardium



ORDER NUMBER	DESCRIPTION	I.D.	O.D.	DILATOR GUIDEWIRE COMPATIBILITY	WORKING LENGTH	OVERALL LENGTH
900005	AcQRef Introducer Sheath	7 F	10.8 F	0.038 in	30 cm	34 cm

For more information please contact your Acutus Medical representative:

ACUTUS MEDICAL, INC.
2210 Faraday Ave
Suite 100
Carlsbad, CA 92008
United States
Phone: +1 442-232-6080
acutus.com

ACUTUS MEDICAL NV
Ikaroslaan 25
1930 Zaventem, Belgium
Phone: +32 2 669 75 00
FAX: +32 2 669 75 01

Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. The AcQMap, AcQGuide and AcQRef devices are CE Marked and FDA cleared.
acutus.com/patents

Acutus Medical®, Acutus®, AcQMap®, AcQGuide® and AcQRef® are registered trademarks of Acutus Medical, Inc. This material is the property of Acutus Medical and its subsidiaries. Please do not copy, forward, or distribute this document.
Copyright © 2020 Acutus Medical, Inc. All rights reserved.

ACUTUS™
M E D I C A L